

GENERAL PURPOSE INCUBATOR INC080-270



GENERAL PURPOSE INCUBATOR INC080-270

INC080-270 GENERAL PURPOSE INCUBATOR

Microprocessor controller (with timing function)

With inner glass door for easy observation

Polished stainless-steel chamber

Independent temperature-limiting alarm system ensures experiments run safely

Printer connector and RS485 connector are options which can connect printer and computer to record the parameters and the variations of temperature



SPECIFICATIONS

Model	INC080-270
Capacity	270 L
Electrical requirement	220 V 50 Hz
Temperature Range	RT+5 - 65°C
Display Resolution	0.1°C/±0.5°C
Ambient Temperature	+5 - 3°C
Temperature Uniformity	±1.5°C (at 37°C)
Power Consumption	750 W
Interior Dimension	600Wx600Dx750H mm
Exterior Dimension	890Wx740Dx910H mm
Shelves	2 pcs
Timing Range	0~5999 min
Remark	With observation window and without inner glass door

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com