

LABSTAC

BIOLOGICAL INDICATOR INCUBATOR INC52-080



labstac.com

BIOLOGICAL INDICATOR INCUBATOR INC52-080

Biological Indicator Incubator is supporting equipment for bio-indicator. It is small in size, less weight and is suitable for pressure steam sterilizer biological indicator, ethylene oxide sterilization biological indicator and other biological indicator sterilization, effect detection and cultivation.

Used in Environmental protection, Research institutes, University laboratories.

Also known as Laboratory Incubator.

INC52-080 BIOLOGICAL INDICATOR INCUBATOR

Designed for incubation of EO biological indicators

Customized wells for placement of biological indicator vials

Easy to set incubation temperature

Aluminium blocks available to meet different experimental requirements



SPECIFICATIONS

Model	INC52-080
Temperature Range	RT+5-80°C
Temperature Accuracy	0.35°C at 57°C
Temperature Uniformity	0.5°C
Temperature Resolution	0.1°C
Time for Heating	≤6 min (20°C to 57°C)
Timer Range	1 min-99h 59 min
Exterior Dimension	110Wx150Dx80H mm
Display	LCD
Weight	0.5/1 kg
Power	120 W
Power Supply	240 V, 60 Hz

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com