

MICROBIOLOGICAL INCUBATOR INC61-070



MICROBIOLOGICAL INCUBATOR INC61-070

Microbiological Incubator can be used for many applications beyond traditional bacterial culture, such as temperature dependent incubations, enzymatic reactions or for reagent storage. Programmable controls and refrigeration are some of the additional features available in incubators on the market today.

Used in Temperature dependent incubations, Enzymatic reactions, Reagent Storage. Programmable controls, Refrigeration.

Also known as Laboratory Microbiological Incubator.

INC61-070 MICROBIOLOGICAL INCUBATOR

Cooling system

High quality imported compressor

Air duct circulation system

Stainless steel inner

Built in BOD socket

large observation window

Lighting device

PID control



SPECIFICATIONS

Model	INC61-070
Capacity	70 L
Temperature Range	0~60°C
Temperature Accuracy	±0.5°C
Temperature Fluctuation	±1%
Timer Range	0-999 min
Internal Dimension	400Wx320Dx500H mm
Exterior Dimension	540Wx570Dx1020H mm
Package Size	620x650x1140 mm
Weight	56.5/62.5 kg
Power	350 W
Power Supply	220 V

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com