

# MICROBIOLOGICAL INCUBATOR INC61-100



# MICROBIOLOGICAL INCUBATOR INC61-100

Microbiological Incubator can be used for many applications beyond traditional bacterial culture, such as temperature dependent incubations, enzymatic reactions or for reagent storage. Programmable controls and refrigeration are some of the additional features available in incubators on the market today.

Used in Temperature dependent incubations, Enzymatic reactions, Reagent Storage. Programmable controls, Refrigeration.

Also known as Laboratory Microbiological Incubator.

## INC61-100 MICROBIOLOGICAL INCUBATOR

Cooling system

High quality imported compressor

Air duct circulation system

Stainless steel inner

Built in BOD socket

Large observation window

Lighting device

PID control



## SPECIFICATIONS

Model	INC61-100
Capacity	100 L
Temperature Range	0~60°C
Temperature Accuracy	±0.5°C
Temperature Fluctuation	±1%
Timer Range	0-999 min
Internal Dimension	450Wx380Dx580H mm
Exterior Dimension	590Wx620Dx1100H mm
Package Size	670x700x1220 mm
Weight	90/98 kg
Power	350 W
Power Supply	220 V

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)