

MICROBIOLOGICAL INCUBATOR INC61-150



MICROBIOLOGICAL INCUBATOR INC61-150

Microbiological Incubator can be used for many applications beyond traditional bacterial culture, such as temperature dependent incubations, enzymatic reactions or for reagent storage. Programmable controls and refrigeration are some of the additional features available in incubators on the market today.

Used in Temperature dependent incubations, Enzymatic reactions, Reagent Storage. Programmable controls, Refrigeration.

Also known as Laboratory Microbiological Incubator.

INC61-150 MICROBIOLOGICAL INCUBATOR

- Cooling system
- High quality imported compressor
- Air duct circulation system
- Stainless steel inner
- Buit in BOD socket
- large observation window
- Lighting device
- PID control



SPECIFICATIONS

Model	INC61-150
Capacity	150 L
Temperature Range	0~60°C
Temperature Accuracy	±0.5°C
Temperature Fluctuation	±1%
Timer Range	0-999 min
Internal Dimension	500Wx380Dx800H mm
Exterior Dimension	640Wx660Dx1350H mm
Package Size	720x740x1470 mm
Weight	102/112 kg
Power	450 W
Power Supply	220 V

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com