

GENERAL PURPOSE INCUBATOR INC81-050



GENERAL PURPOSE INCUBATOR INC81-050

General Purpose Incubator are designed to accommodate an ample range of laboratory applications and lends a contamination-free environment for cell and tissue cultures. These incubators serve the purpose for most biological analysis, research, general laboratory applications, biochemical studies, bacteriology, and hematological studies. Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing. Also known as Lab Incubator, Laboratory General Purpose Incubator.

INC81-050 GENERAL PURPOSE INCUBATOR

LCD PID control system

International brand cooling fan

Stainless steel chamber

Cavity preheating technology and duct cycling system



SPECIFICATIONS

Model	INC81-050
Capacity	50 L
Temperature Range	RT +5 ~ 65°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±0.5%
Internal Dimension	380Wx380Dx350H mm
Exterior Dimension	580Wx530Dx600H mm
Shelves	2
Weight	31 kg
Power Supply	220 V 50 Hz

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com