

# GENERAL PURPOSE INCUBATOR INC82-150



# GENERAL PURPOSE INCUBATOR INC82-150

General Purpose Incubator are designed to accommodate an ample range of laboratory applications and lends a contamination-free environment for cell and tissue cultures. These incubators serve the purpose for most biological analysis, research, general laboratory applications, biochemical studies, bacteriology, and hematological studies. Used in Quality Control, Medicine Processing Industries, Research, Pharmaceuticals, Microbial Culturing. Also known as Lab Incubator, Laboratory General Purpose Incubator.

## INC82-150 GENERAL PURPOSE INCUBATOR

- Balance Refrigeration
- Adopt without fluorine R134a refrigerant technology
- Duct circulation
- Stainless steel inner
- Programmable intelligent PID control
- Usd large capacity humidifier
- Double door design
- Lighting device
- With UV lamp



### SPECIFICATIONS

Model	INC82-150
Capacity	150 L
Temperature Range	0-60°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±1%
Humidity Range	50%-90 %
Humidity Fluctuation	±5%
Timer Range	0-999 min
Internal Dimension	500Wx380Dx800H mm
Exterior Dimension	640Wx660Dx1350H mm
Package Size	720x740x1470 mm
Weight	105/115 kg
Power	600 W
Power Supply	220 V

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)