

ANAEROBIC INCUBATOR INCD1-362



ANAEROBIC INCUBATOR INCD1-362

Anaerobic Incubator is a workstation designed for various studies of culture anaerobic micro-organisms. Used in Used to culture anaerobic micro-organism to perform various studies in pathology labs, soil research, oil industry.

Also known as Anaerobic Incubator.

INCD1-362 ANAEROBIC INCUBATOR

Equipped with large LCD screen display, high-precision microcomputer control (with Timer function), it can accurately and visually reflect the actual temperature inside the incubator, over heating alarm adopted, safe and reliable.

Switch control for solenoid valve, it can accurately adjust the flow and input any necessary gas.

UV Sterilizer, effectively prevent bacterial contamination.

Stainless steel cultivation and operation room, transparent impact-resistant glass front window for easy observation. Latex gloves are comfortable and reliable, easy to use.

Operation room is equipped with deoxidization catalyst.

Equipped with leakage protection.

The incubator is designed with a double wide door, it can put more petri dishes.



Product Image Coming Soon

SPECIFICATIONS

Model	INCD1-362
Capacity	< 5 minutes
Oxygen	362 L
RCF. Range	< $\pm 0.3^{\circ}\text{C}$
Gross Capacity	> 12hrs (when no supply of mixed gas)
Temperature Range	RT+3~60°C

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com