

AUTOMATIC SEALING AND CAPPING MACHINE LIH11-16S



AUTOMATIC SEALING AND CAPPING MACHINE

LIH11-16S

Automatic Sealing and Capping Machine is specifically designed to speed up the process of packaging of your product but also ensure the sterility of the product and that the product does not get contaminated

Used in Pharmaceuticals, Mobile inspection laboratories, Disease prevention and control centers, Food Industry, Beverage Industry.

Also known as Screw on capping machine.

LIH11-16S AUTOMATIC SEALING AND CAPPING MACHINE


Fully automatic and one key operation, no need of additional manual operation

Adapts to varieties of plate types

Adjustable sealing pressure

A variety of communication interfaces are reserved, which can be quickly connected to other machines

Each roll of film can seal more than 350 plates(reference quantity)



SPECIFICATIONS

Model	LIH11-16S
Positions	1pc
Method	Cold Seal
Speed	<16s/plate
Operating Mode	Manual plate change, automatic sealing
Application	PCR plate, deep well plate and micro-plate conforming to SBS standard(Customizable)
Control Mode	Buttons
Interface	USB
External Size	310x380x410 mm
Package Size	454x382x464 mm
Weight	Net Weight- 16.5kg Gross Weight- 20.5kg
75 W	Power Supply

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com