

**LABSTAC**

# **AUTOMATIC SEALING AND CAPPING MACHINE LIH11-96T**



# AUTOMATIC SEALING AND CAPPING MACHINE

## LIH11-96T

Automatic Sealing and Capping Machine is specifically designed to speed up the process of packaging of your product but also ensure the sterility of the product and that the product does not get contaminated  
Used in Pharmaceuticals, Mobile inspection laboratories, Disease prevention and control centers, Food Industry, Beverage Industry.  
Also known as Screw on capping machine.

### LIH11-96T AUTOMATIC SEALING AND CAPPING MACHINE

- One-button start
- Simple and fast operation
- Safe and labor-saving
- Capping/sealing speed  $\leq 4s$



#### SPECIFICATIONS

Model	LIH11-96T
Throughput	8-96, suitable for 8-strip PCR tubes and 96-well PCR plates
Start Method	Button
Speed	Capping\sealing speed $\leq 4s$
Function	8-strip PCR tubes capping, 96-well PCR plate capping, PCR plate sealing film
Adapter	2 types, 8-strip PCR tubes and 96-well PCR plates, can be customized
Pressure	Peak pressure: 300kg
External Size	277x194x230 mm
Package Size	300x260x340 mm
Weight	Net Weight- 9.2kg Gross Weight- 12kg
Power Supply	100-240V 50/60Hz



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)