

# PORTABLE ULTRASONIC PROCESSOR LIH14-100



# PORTABLE ULTRASONIC PROCESSOR LIH14-100

Portable Ultrasonic Processor is used to emulsify, disperse, homogenization, extraction, anti-foaming, cleaning, preparation of nano-materials, dispersion, and speed up chemical reactions

Used in Biotechnology, Pharmaceutical, Chemical, Industrial, Environmental.

Also known as Portable Ultrasonic Homogenizer Sonicator, Portable ultrasonic cell crusher.

## LIH14-100 PORTABLE ULTRASONIC PROCESSOR



### SPECIFICATIONS

Model	LIH14-100
Operating frequency range	20-25 KHz, frequency automatic tracking
The end of the horn diameter	Φ 3 or 5
Ultrasonic Power	150W, can be adjusted, hit into ways to fine-tune each of 5W
Parameters sample handling capacity	0.1-100ml (need to select the appropriate probe)
Generator output power	150W (max)
Operating voltage	220V±5%

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)