

PORTABLE ULTRASONIC PROCESSOR LIH14-100A



PORTABLE ULTRASONIC PROCESSOR LIH14-100A

Portable Ultrasonic Processor is used to emulsify, disperse, homogenization, extraction, anti-foaming, cleaning, preparation of nano-materials, dispersion, and speed up chemical reactions
Used in Biotechnology, Pharmaceutical, Chemical, Industrial, Environmental.
Also known as Portable Ultrasonic Homogenizer Sonicator, Portable ultrasonic cell crusher.

LIH14-100A PORTABLE ULTRASONIC PROCESSOR



SPECIFICATIONS

Model	LIH14-100A
Operating frequency range	20-25 KHz, frequency automatic tracking
The end of the horn diameter	Φ 3 or 5
Ultrasonic Power	250W, can be adjusted, hit into ways to fine-tune each of 5W
Parameters sample handling capacity	0.1-100ml (need to select the appropriate probe)
Generator output power	250W (max)
Operating voltage	220V±5%

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com