

MULTI-PARAMETER ANALYZER MET101-074



MULTI-PARAMETER ANALYZER MET101-074

Multi-parameter Analyzer simultaneously measures Ion, Conductivity, Dissolved Oxygen, Temperature that otherwise requires four different meters. With the appropriate probes installed, it can also simultaneously measure several or one parameters.

MET101-074 MULTI-PARAMETER ANALYZER



Product Image Coming Soon

SPECIFICATIONS

Model	MET101-074
DC Voltage	
200 mV/2 V/20 V/200 V	± (0.5%+2)
600 V	± (0.8%+2)
AC Voltage	
2 V/20 V/200 V	± (1.0%+3)
600 V	± (1.2%+3)
DC Current	
200 µA/2 mA/20 mA/200 mA	± (1.5%+3)
AC Current	
200 µA/2 mA/20 mA/200 mA	± (1.5%+4)
DC	
10 A	± (3.0%+3)
AC	
10A	± (3.0%+3)
Resistance	
200 Ω/2 kΩ/20 kΩ/200 kΩ/2 MΩ	± (0.8%+2)
20MΩ	± (1.0 % + 2)
Temperature	
- 20 °C ~ 1000 °C	± (3.0 % + 3)
0 °F ~ 1800 °F	± (3.0 % + 3)
Display	2000 Counts
Data Hold	Yes
Continuity	Yes
Diode Test	Yes

Maximum Value Hold	Yes
Power Supply	1x9 V 6F22 Battery
Weight	Approx. 172.5/173 g
Size	148x74x42 mm
Safety Standard/Rating	ENG1010-1,ENG1010-2-033, EN61326,CAT III 600 V

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com