

# VERTICAL ROTATING MIXER MIX12-80V



# VERTICAL ROTATING MIXER MIX12-80V

DC brushless motor provides reliable and quiet running, maintenance free and can be used in cold rooms or in incubators, Majorly used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology and Wastewater Facilities.

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Tube Rotator, Rotator, Laboratory Rotator.

## MIX12-80V VERTICAL ROTATING MIXER

Smart appearance, solid structure, simple operation

LED display time and speed

DC brushless motor, stable performance, low noise and long using life

360° vertical circulation rotation, operating smooth and steady

Different tube holder frames for optional to suit different types of tubes



## SPECIFICATIONS

Model	MIX12-80V
Speed Range	10 ~ 80 rpm
Shaking Orbit	circle
Time Range	1min~99h59min
Voltage	AC 100~230V, 50/60Hz
Power	16 W
Fuse	250V, 1A, $\Phi$ 5x20
Dimension	W.470 x D.190 x H.160mm
Net Weight	2.8 kgs

## ACCESSORIES

Accessory Code	Name	Capacity
3200906005	Tube Holder frame	1.5mlx64
3200906006	Tube Holder frame	15mlx22
3200906007	Tube Holder frame	50mlx16

***LABSTAC***

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)