

ROTARY MIXING SOLUTION MIX13-20



ROTARY MIXING SOLUTION MIX13-20

Constructed of rugged, corrosion-resistant materials, the instruments endure heavy everyday usage. The ergonomic low profile design makes everyday vortexing comfortable for the user. Unnecessary movement has been avoided during use due to robust die-cast body.

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

MIX13-20 ROTARY MIXING SOLUTION

It can be used in the humid environment such as refrigerator and refrigerated cabinet

Several optional turntable supplied, suit for the tubes from 1.5ml to 50ml

One button to start/stop, easy to operate, can run for a long time, running smoothly without noise at 20 rpm



SPECIFICATIONS

Model	MIX13-20
Speed Range	20 rpm
Gradient Angle Range	45° (fixed angle)
Voltage	AC 220V/AC 110V, 50/60 Hz
Power	14 W
Fuse	250V, 1A, Φ 5x20
Dimension	W.260 x D.300 x H.233mm
Net Weight	1.8 kgs

OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Capacity
3200907005	Turntable for rotary mixing solution	50ml V-bottom tube x 10	
3200907006	Turntable for rotary mixing solution	1.5ml tube x 64	
3200907007	Turntable for rotary mixing solution	15ml V-bottom tube x 18	
3200907008	Foam tube holder		50ml centrifuge tubex12 15ml centrifuge tubex12 1.5ml centrifuge tubex30
3200907009	Turntable		50ml centrifuge tubex10 15ml centrifuge tubex10 1.5ml centrifuge tubex30
3200907010	Double turntable connector		

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com