

**LABSTAC**

# VERTICAL ROTATING MIXER MIX13-80G



labstac.com

# VERTICAL ROTATING MIXER MIX13-80G

DC brushless motor provides reliable and quiet running, maintenance free and can be used in cold rooms or in incubators, Majorly used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology and Wastewater Facilities.

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Tube Rotator, Rotator, Laboratory Rotator.

## MIX13-80G VERTICAL ROTATING MIXER

Humanized designed operation keyboard

LED display time and speed

Several optional turntable supplied, suit for the tubes from 1.5ml to 50ml

The turntable's gradient can be adjusted from 0 to 90 degree

Timing function

Beep alarm when program time ends

DC brushless motor, stable performance, low noise and long using life



## SPECIFICATIONS

Model	MIX13-80G
Speed Range	10 ~ 80 rpm
Gradient Angle Range	0°~90°
Time Range	1min~9999min
Voltage	AC 100~230V, 50/60 Hz
Power	35 W
Fuse	250V, 1A, $\Phi$ 5x20
Dimension	W.206 x D.260 x H.300mm
Net Weight	4.5 kgs

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description	Capacity
3200907005	Turntable for rotary mixing solution	50ml V-bottom tube x 10	
3200907006	Turntable for rotary mixing solution	1.5ml tube x 64	
3200907007	Turntable for rotary mixing solution	15ml V-bottom tube x 18	
3200907008	Foam tube holder		50ml centrifuge tubex12 15ml centrifuge tubex12 1.5ml centrifuge tubex30
3200907009	Turntable		50ml centrifuge tubex10 15ml centrifuge tubex10 1.5ml centrifuge tubex30
3200907010	Double turntable connector		

***LABSTAC***

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)