

TUBE ROLLER MIX14-80



TUBE ROLLER MIX14-80

DC brushless motor provides reliable and quiet running, maintenance free and can be used in cold rooms or in incubators, Majorly used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology and Wastewater Facilities.

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Tube Rotator, Rotator, Laboratory Rotator.

MIX14-80 TUBE ROLLER

End-over-end or rolling action, gives gentle but effective mixing

Adjustable speed range of 0 to 80 rpm

Used with a wide variety of accessories

Continuous quiet operation



SPECIFICATIONS

Model	MIX14-80
Speed Range	0-80 rpm
Permissible Relative Humidity	0.8
Tilt Angle (°)	90
Operation	Continuous
Dimension	150Dx530Wx190H mm
Motor Type	DC motor
Protection Class	IP21
Weight	1.7 kg
Power	20 W
Power Supply	100-240V, 50/60Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Description
3200705005	Disk and Rotisserie accessory	For 1.5ml x 12 centrifuge tubes
3200705006	Disk and Rotisserie accessory	For 15ml x 12 centrifuge tubes
3200705007	Disk and Rotisserie accessory	For 50ml x 12 centrifuge tubes

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com