

# TUBE ROLLER MIX15-40



# TUBE ROLLER MIX15-40

DC brushless motor provides reliable and quiet running, maintenance free and can be used in cold rooms or in incubators, Majorly used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology and Wastewater Facilities.

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Tube Rotator, Rotator, Laboratory Rotator.

## MIX15-40 TUBE ROLLER

The rotator comes with a disk type design

Digital control speed and allowing procedures to be accurate

LCD displays timer


End-over-end or rolling action, gives gentle but effective mixing

Adjustable speed range of 10 to 70rpm

Adjustable mixing angle between 0 and 90 °

Used with a wide variety of accessories

Continuous and timed quiet operation



## SPECIFICATIONS

Model	MIX15-40
Speed Range	10-70 rpm
Permissible Ambient Temperature	5-40°C
Permissible Relative Humidity	0.8
Tilt Angle (°)	0-90
Operation	Continuous/ timed operation
Timer	Yes
Timing Range	1-1199 min
Dimension	510Dx190Wx260H mm
Motor Type	DC motor
Protection Class	IP21
Display	LCD
Weight	8 kg
Power	40 W
Power Supply	100-240V, 50/60Hz

## OPTIONAL ACCESSORIES

Accessory Code	Name	Description
3200706005	Disk accessory	For 1.5ml x 60 centrifuge tubes holder
3200706006	Disk accessory	For 15ml x 16 centrifuge tubes holder
3200706007	Disk accessory	For 50ml x 8 centrifuge tubes holder

3200706008	Disk support rods	4 pcs; use with BX-103920; BX-103921; BX-103922
------------	-------------------	---



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)