

VORTEX MIXER MIX31-28C



VORTEX MIXER MIX31-28C

Constructed of rugged, corrosion-resistant materials, the instruments endure heavy everyday usage. The ergonomic low profile design makes everyday vortexing comfortable for the user. Unnecessary movement has been avoided during use due to robust die-cast body.

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

MIX31-28C VORTEX MIXER

The perfect mixing radius combined with excellent 2D mixing technology and adjustable mixing speed ensures fast and perfect mixing of multiple samples.



SPECIFICATIONS

Model	MIX31-28C
Speed Range	300 ~3000 rpm
Touch Vortex Speed	3500 rpm
Shaking Orbit	3 mm
Time Range	1 s~99 h 59 min
Working Noise	<50dB
Voltage	AC 100~230V, 50/60Hz
Power	45 W
Fuse	250 V, 1 A, φ5 x 20
Dimension	W.196 x D.270 x H.170 mm
Net Weight	7.7 kgs

ACCESSORIES

Accessory Code	Name	Capacity
3000510005	Tube Stand	96x0.2ml tube stand. Applied to 0.2ml PCR tubes; PCR strips and 96-wellPCR plates(semi-skirted and unskirted)
3000510006	Tube Stand	4x0.5ml Tube stand
3000510007	Tube Stand	24x1.5/2.0ml Tube stand

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com