

# MICROPLATE READER MRW11-01



# MICROPLATE READER MRW11-01

Microplate Reader is a laboratory instrument that is used to measure chemical, biological or physical reactions, properties and analytes within the well of a microplate. Flexibility and reliability makes this product an unique choice. Used in Life Science, Drug Discovery, Research, Laboratory, Medical, Assay Development.  
Also known as Laboratory Microplate Reader.

## MRW11-01 MICROPLATE READER

It is an eight-channel reader with the most sophisticated technology, high quality, high reliability and easy to operate.

8 channels measuring system.

Analysis Mode: Cut-off, single and multi-standard, O.D., bichromatic, dynamic reading.

Optical filter 405/450/492/630nm

100 programs storage.

Touch screen display.

Built-in thermal printer.



## SPECIFICATIONS

Model	MRW11-01
Measuring System	Optical system with 8 channels
Analysis Mode	Cut-off, single and multi-standard, O.D., bichromatic, dynamic reading
Linear Range	0.001-3.500 Abs
Photometric Accuracy	±1% or ±0.001 Abs
Repeatability	<1.0%
Stability	<0.005 Abs/hr
Linearity	<±1.0%
Sensitivity	≤±0.001A
Channels Discrepancy	≤±0.002A
Optical Filters	405/450/492/630 nm other filters available on request
Reading	Automatic focusing function makes the results more precise
Tray Shaking	Yes
Storage	100 programs
Display	Touch screen
Interface	RS-232 and USB port
Printer	Built-in thermal printer (external printer are available on request)
Power Supply	AC 220V 50 Hz or 110V 60 Hz
Dimension	46Lx33Wx18H cm
Weight	7.5 kg

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)