

MICROCIRCULATION MICROSCOPE MSC11-0480



MICROCIRCULATION MICROSCOPE MSC11-0480

Micro-circulation Microscope is an advanced medical Instrument equipped with a special cold light source, all kinds of microscopic examination of the capillaries of the skin, mucosa and muscle are trauma-less. It can be used for ecological observation, photographing and video-recording the micro-circulatory status.

Used in Microbiology Lab, Pharmaceutical, Hospitals, Medical, Universities, Scientific research, Institute.

Also known as Ultra microscope, High-magnification microscope, High resolution microscope.

MSC11-0480 MICROCIRCULATION MICROSCOPE

Excellent optical function with infinite optical system.


Innovative stand structure and sharp image display.

Valid mould proof design, convenient operation and special design for viewing incubating cell tissue.

It can be used under the high temperature and high humidity environment.

Focus System: Coaxial coarse/ Slightly adjusted, Z axis movement, Mechanical upper and lower limits.

Illumination: LED illumination with high brightness and lower consumption, new luminaire with lower temperature and 10,000-200,000h working life, 90°revolve left and right.



Product Image Coming Soon

SPECIFICATIONS

Model	MSC11-0480
Magnification	480X
Display Size	204x153 mm
Objective	10X
Stage	360°revolve, X,Y axis movement ,moving range 30x80 mm
Camera	Color camera system, HD LCD display
Package Size	300x380x380 mm
Gross Weight	7 kg
Power Supply	12 V

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com