

DOUBLE ARM LED OPERATING LAMP OPL15-14



DOUBLE ARM LED OPERATING LAMP OPL15-14

Challenging tasks require special solutions. The innovative Labstac LED operation light give you benefits of cutting-edge LED technology in a uniquely effective manner. We focus on providing LED lights that are more color-accurate, more flexible, more economical. We offer innovative light and system solutions that can be retrofitted easily and without special preparation.

Used in Dental implant surgery, Medical unit, Surgery, Operation theatre.

Also known as Surgical LED Light, LED operation theatre light, LED OT light.

OPL15-14 DOUBLE ARM LED OPERATING LAMP

Streamlined, ultra-thin, enclosed head with reduced air resistance to minimise contamination

Digital controlled light intensity without grades of brightness and automatic memory for added flexibility

Optionally built in power supply systems avoiding complicated installation of wall mounted control boxes

Auxiliary light source included for laparoscopic operations

Rigid suspension arm with easy mobility and stable positioning


360° arm rotation

No infrared or ultraviolet radiation

Backlight compensation function

Endoscopy light source

ONDAL arm, OSRAM bulb



SPECIFICATIONS

Model	OPL15-14
Illumination	140,000 Lux
Adjusting range of illuminance	10-100%
Spot diameter (D50-D10)	130-350 mm
Color rendering Index	≥97
Color temperature	3600K-5500K (±500K) variable
Average bulb life (hours)	100,000 h Power of bulb- 1 W x N 64 Light focus system- Manual
Power Supply	220 V 50 Hz
Input power	80 W
Ceiling height for installation	2.8 m

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com