

# HOLE-TYPE OR LAMP OPL39-50I



# HOLE-TYPE OR LAMP OPL39-50I

Hole type operation light is suitable for small operating room. Labstac lights uses multiple different sources and by focusing knob, the beam is illuminated from different surgery locations which reduces the shadow by medical staff head, hands and other parts caused. The illumination adjusted according to actual needs.

Used in Medical unit, Surgery, Operation theatre.

Also known as Hole type surgical Light, Hole type OT light.

## OPL39-50I HOLE-TYPE OR LAMP

Combined type lighthead



### SPECIFICATIONS

Model	OPL39-50I
Illumination	80,000LUX (at 1m distance)
Color temperature	4500K±500K
Color rendering Index	100≥Ra≥85
Depth of light beam	>500 cm
The Value of total irradiance	<6 mW/m <sup>2</sup> lx
Total irradiance	<600 W/m <sup>2</sup>
Temperature rise on the head of the patient under surgical operation	≤2°C
Rated Voltage of Bulb	24V
Rated Power of bulb	125 W
Input power	150 VA
Power Supply	220V±10% 50Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)