

HOLE-TYPE OR LAMP OPL3B-12



HOLE-TYPE OR LAMP OPL3B-12

Hole type operation light is suitable for small operating room. Labstac lights uses multiple different sources and by focusing knob, the beam is illuminated from different surgery locations which reduces the shadow by medical staff head, hands and other parts caused. The illumination adjusted according to actual needs.

Used in Medical unit, Surgery, Operation theatre.

Also known as Hole type surgical Light, Hole type OT light.

OPL3B-12 HOLE-TYPE OR LAMP



SPECIFICATIONS

| | |
|--|-----------------------------------|
| Model | OPL3B-12 |
| Illumination | $\geq 160,000\text{LUX}$ |
| Color temperature | $4500\text{K}\pm 500\text{K}$ |
| Color rendering Index | $100\geq Ra\geq 85$ |
| Depth of light beam | $>500\text{ cm}$ |
| The Value of total irradiance | $<6\text{ mW/m}^2\text{lx}$ |
| Total irradiance | $<600\text{ W/m}^2$ |
| Temperature rise on the head of the patient under surgical operation | $\leq 2^{\circ}\text{C}$ |
| Rated Voltage of Bulb | 24V |
| Rated Power of bulb | 300 W |
| Input power | 350 VA |
| Power Supply | $220\text{V}\pm 10\% 50\text{Hz}$ |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com