

MANUAL OR TABLE OPT12-2000B



MANUAL OR TABLE OPT12-2000B

Manual Operation table are suitable for various operation tables conducted by medical unit. The OR table and accessories are ergonomic and can be used across multiple disciplines. Various actions are controlled by hydraulic pneumatic system. 360° rotation helps every surgical requirement.

Used in Medical unit, Surgery, Operation theatre, cardiovascular, pediatric, gynecology, abdominal, gall bladder, plastic surgery and more.

Also known as Manual surgical table, Manual OR Table, Manual operating table.

OPT12-2000B MANUAL OR TABLE

Imported LINAK motor various actions are controlled by remote controller

The table, base seat and upright column covering are made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned

The mattress adopts high density memory sponge and is once formed, with no seam

Whole stainless steel hidden auxiliary table with flexible extension and retraction

Equipped with stainless steel sewage basin and other accessories

Foot treadle brake device which is easy for moving obstetric table



SPECIFICATIONS

Model	OPT12-2000B
Table Length	2000±50 mm
Table Width	600±20 mm
Height Adjustment	(700~950) ±50 mm
Working area length/width	520x500 mm
Trendelenburg/Rev-Trend	≥10°/≥12°
Backboard upwards	≥50°
Standard Accessory	A pair of arm support, a pair of leg support, a pair of foot pedal, a pair of arm rods, a stainless steel sewage basin, a handheld remote, controller
Main Voltage AC	220V±22V, 50Hz±1Hz
Input Power VA	150±10%

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com