

# ORTHOPEDICS TRACTOR RACK OPT14-1700



# ORTHOPEDICS TRACTOR RACK OPT14-1700

Labstac's Orthopedics tractor rack is mainly applied to orthopedics surgical by using this tractor rack the patient's leg can be performed side traction, prone traction and open-up position of the body during the operation.

Used in Surgery, Leg operation, Medical unit.

Also known as Operating Table.

## OPT14-1700 ORTHOPEDICS TRACTOR RACK

Adopt carbon fiber material, the prespective imaging is clear, light and strong

The plug in design is convient and quick to connect with table

Equipped with a traction frame accessory car, which is convient to move and save space



## SPECIFICATIONS

Model	OPT14-1700
Table Length	1700±10 mm
Table Width	750±100 mm
Traction device horizontal rotary	0-180°
Traction distance	(0-200)+ 10%

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)