

DRYING OVEN OVE075-760



DRYING OVEN OVE075-760

OVE075-760 DRYING OVEN

304 stainless steel, mirror polishing processing, semicircular arcs at corners easy to clean and maintain. The space between the shelves in the chamber is adjustable

PID controller with over temperature alarm and timing function ensures precise and reliable control

Ceramic fiber door seal, which can run at high temperature for a long time and has a long service life

Programmable controller: 7 periods 63 steps, 0 ~5999mins for each period, can preset boot and shutdown time, adjustable circulating fan

Independent over-temperature alarm system ensures experiments running safely



SPECIFICATIONS

Model	OVE075-760
Capacity	760 L
Electrical requirement	AC 380 V 50 HZ
Temperature Range	RT+20°C 400°C
Display Resolution	0.1°C
Temperature Stability	±0.5°C
Ambient Temperature	+5 - 40°C
Power Consumption	12750 W
Internal Dimension	980Wx1000Dx780H mm
External Dimension	1324Wx1263Dx1770H mm
Shelves	2 pcs
Timing Range	0~5999 min

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com