

VACUUM DRYING OVEN OVE076-024C



VACUUM DRYING OVEN OVE076-024C

OVE076-024C VACUUM DRYING OVEN

304 stainless steel, mirror polishing processing, easy to clean and maintain
PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time
A big dual layer tempered glass on the door provides good observation
Door adjustment system with silicon door gasket ensures better vacuity



SPECIFICATIONS

Model	OVE076-024C
Capacity	24 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	700W
Temperature range	RT-10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	300Wx300Dx275H mm
Exterior Dimension	605Wx490Dx450H mm
Shelves	4 pcs
Chamber Material	304 Stainless steel

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com