


# VACUUM DRYING OVEN OVE076-030A



# VACUUM DRYING OVEN OVE076-030A

## OVE076-030A VACUUM DRYING OVEN

304 stainless steel, mirror polishing processing, easy to clean and maintain  
PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time  
A big dual layer tempered glass on the door provides good observation  
Door adjustment system with silicon door gasket ensures better vacuity



### SPECIFICATIONS

Model	OVE076-030A
Capacity	30 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	550 W
Temperature range	RT-10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	320Wx320Dx300H mm
Exterior Dimension	630Wx510Dx490H mm
Shelves	1 pcs
Chamber Material	316 Stainless steel



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)