

VACUUM DRYING OVEN OVE076-030C



VACUUM DRYING OVEN OVE076-030C

OVE076-030C VACUUM DRYING OVEN

304 stainless steel, mirror polishing processing, easy to clean and maintain
PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time
A big dual layer tempered glass on the door provides good observation
Door adjustment system with silicon door gasket ensures better vacuity



SPECIFICATIONS

Model	OVE076-030C
Capacity	30 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	850 W
Temperature range	RT-10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	320Wx320Dx300H mm
Exterior Dimension	630Wx510Dx490H mm
Shelves	4 pcs
Chamber Material	316 Stainless steel

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com