

VACUUM DRYING OVEN OVE076-050D



VACUUM DRYING OVEN OVE076-050D

OVE076-050D VACUUM DRYING OVEN

304 stainless steel, mirror polishing processing, easy to clean and maintain
PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time
A big dual layer tempered glass on the door provides good observation
Door adjustment system with silicon door gasket ensures better vacuity



SPECIFICATIONS

Model	OVE076-050D
Capacity	50 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	1450 W
Temperature range	RT-10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	415Wx370Dx345H mm
Exterior Dimension	730Wx560Dx550H mm
Shelves	6 pcs
Chamber Material	304 Stainless steel

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com