

# VACUUM DRYING OVEN OVE076-050E



# VACUUM DRYING OVEN OVE076-050E


## OVE076-050E VACUUM DRYING OVEN

304 stainless steel, mirror polishing processing, easy to clean and maintain

PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time

A big dual layer tempered glass on the door provides good observation

Door adjustment system with silicon door gasket ensures better vacuity



## SPECIFICATIONS

Model	OVE076-050E
Capacity	50 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	650 W
Temperature range	RT-10 - 65°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	415Wx370Dx345H mm
Exterior Dimension	730Wx560Dx550H mm
Shelves	2 pcs
Chamber Material	304 Stainless steel

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)