

# VACUUM DRYING OVEN OVE076-090B



# VACUUM DRYING OVEN OVE076-090B

## OVE076-090B VACUUM DRYING OVEN

304 stainless steel, mirror polishing processing, easy to clean and maintain  
PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time  
A big dual layer tempered glass on the door provides good observation  
Door adjustment system with silicon door gasket ensures better vacuity



### SPECIFICATIONS

Model	OVE076-090B
Capacity	90 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	870 W
Temperature range	RT-10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	450Wx450Dx450H mm
Exterior Dimension	740Wx610Dx591H mm
Shelves	4 pcs
Chamber Material	304 Stainless steel

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)