

# VACUUM DRYING OVEN OVE076-125B



# VACUUM DRYING OVEN OVE076-125B

## OVE076-125B VACUUM DRYING OVEN

304 stainless steel, mirror polishing processing, easy to clean and maintain  
PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time  
A big dual layer tempered glass on the door provides good observation  
Door adjustment system with silicon door gasket ensures better vacuity



### SPECIFICATIONS

Model	OVE076-125B
Capacity	125 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	1500 W
Temperature range	RT-10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	500Wx500Dx500H mm
Exterior Dimension	790Wx660Dx641H mm
Shelves	6 pcs
Chamber Material	304 Stainless steel



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)