

VACUUM DRYING OVEN OVE076-216A



VACUUM DRYING OVEN OVE076-216A


OVE076-216A VACUUM DRYING OVEN

304 stainless steel, mirror polishing processing, easy to clean and maintain

PID controller with over temperature alarm and timing function ensures precise and reliable control, also save more than 40% heating time

A big dual layer tempered glass on the door provides good observation

Door adjustment system with silicon door gasket ensures better vacuity



SPECIFICATIONS

Model	OVE076-216A
Capacity	216 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	1550 W
Temperature range	RT-10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	600Wx600Dx600H mm
Exterior Dimension	890Wx760Dx741H mm
Shelves	3 pcs
Chamber Material	304 Stainless steel

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com