

VACUUM DRYING OVEN OVE077-090B



VACUUM DRYING OVEN OVE077-090B

OVE077-090B VACUUM DRYING OVEN

- Microprocessor controller
- Dual layer tempered glass door for clear observation
- Minimum heating time 50% less than traditional vacuum oven
- Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning
- Standard with vacuum pump



SPECIFICATIONS

Model	OVE077-090B
Capacity	90 L
Electrical requirement	AC 220 V 50 HZ
Controller	LED display
Power Consumption	1350 W
Temperature range	RT+10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	450Wx450Dx450H mm
Exterior Dimension	610Wx590Dx1350H mm
Shelves	4 pcs
Chamber Material	304 Stainless steel



Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com