

# VACUUM DRYING OVEN OVE077-215B



# VACUUM DRYING OVEN OVE077-215B

## OVE077-215B VACUUM DRYING OVEN

Microprocessor controller

Dual layer tempered glass door for clear observation

Minimum heating time 50% less than traditional vacuum oven

Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning

Standard with vacuum pump



## SPECIFICATIONS

|                        |                     |
|------------------------|---------------------|
| Model                  | OVE077-215B         |
| Capacity               | 215 L               |
| Electrical requirement | AC 220 V 50 HZ      |
| Controller             | LED display         |
| Power Consumption      | 2100 W              |
| Temperature range      | RT+10 - 200°C       |
| Display Resolution     | 0.1                 |
| Temperature Stability  | ±1°C                |
| Vacuum Degree          | 133 Pa              |
| Ambient Temperature    | +5 - 40°C           |
| Interior Dimension     | 560Wx600Dx640H mm   |
| Exterior Dimension     | 720Wx8200Dx1750H mm |
| Shelves                | 6 pcs               |
| Chamber Material       | 304 Stainless steel |

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)