

VACUUM DRYING OVEN OVE077-431



VACUUM DRYING OVEN OVE077-431

OVE077-431 VACUUM DRYING OVEN

Microprocessor controller

Dual layer tempered glass door for clear observation

Minimum heating time 50% less than traditional vacuum oven

Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning

Standard with vacuum pump

Direct heating and control by shelves: Each shelf has a separate temperature sensor, can be independently displayed and controlled temperature in setting, monitoring and timing, more accurate and stable



SPECIFICATIONS

| | |
|------------------------|----------------------|
| Model | OVE077-431 |
| Capacity | 431 L |
| Electrical requirement | AC 380 V 50 HZ |
| Controller | LED display |
| Power Consumption | 3800 W |
| Temperature range | RT+10 - 200°C |
| Display Resolution | 0.1 |
| Temperature Stability | ±1°C |
| Vacuum Degree | 133 Pa |
| Ambient Temperature | +5 - 40°C |
| Interior Dimension | 630Wx810Dx845H mm |
| Exterior Dimension | 1000Wx1040Dx1855H mm |
| Shelves | 4 pcs |
| Chamber Material | 304 Stainless steel |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com