

VACUUM DRYING OVEN OVE077-431



VACUUM DRYING OVEN OVE077-431

OVE077-431 VACUUM DRYING OVEN

- Microprocessor controller
- Dual layer tempered glass door for clear observation
- Minimum heating time 50% less than traditional vacuum oven
- Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning
- Standard with vacuum pump
- Direct heating and control by shelves: Each shelf has a separate temperature sensor, can be independently displayed and controlled temperature in setting, monitoring and timing, more accurate and stable



SPECIFICATIONS

Model	OVE077-431
Capacity	431 L
Electrical requirement	AC 380 V 50 HZ
Controller	LED display
Power Consumption	3800 W
Temperature range	RT+10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	630Wx810Dx845H mm
Exterior Dimension	1000Wx1040Dx1855H mm
Shelves	4 pcs
Chamber Material	304 Stainless steel



Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com