

VACUUM DRYING OVEN OVE078-125



VACUUM DRYING OVEN OVE078-125

OVE078-125 VACUUM DRYING OVEN

- Microprocessor controller
- Dual layer tempered glass door for clear observation
- Minimum heating time 50% less than traditional vacuum oven
- Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning
- Standard with vacuum pump
- Provides programmability for 7 periods and 9 steps for each period making for a total of 63 programmable step



SPECIFICATIONS

Model	OVE078-125
Capacity	125 L
Electrical requirement	AC 220 V 50 HZ
Controller	LCD Programmability controller
Power Consumption	2050 W
Temperature range	RT+10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	500Wx500Dx500H mm
Exterior Dimension	660Wx640Dx1400H mm
Shelves	3 pcs
Chamber Material	304 Stainless steel

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com