

# VACUUM DRYING OVEN OVE078-913



# VACUUM DRYING OVEN OVE078-913

## OVE078-913 VACUUM DRYING OVEN

Microprocessor controller

Dual layer tempered glass door for clear observation

Minimum heating time 50% less than traditional vacuum oven

Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning

Standard with vacuum pump

Provides programmability for 7 periods and 9 steps for each period making for a total of 63 programmable step



## SPECIFICATIONS

Model	OVE078-913
Capacity	913 L
Electrical requirement	AC3 80 V 50 HZ
Controller	LCD Programmability controller
Power Consumption	5600 W
Temperature range	RT+10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	750Wx1160Dx1050H mm
Exterior Dimension	1400Wx1395Dx2010H mm
Shelves	5 pcs
Chamber Material	304 Stainless steel

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)