

VACUUM DRYING OVEN OVE078-913



VACUUM DRYING OVEN OVE078-913

OVE078-913 VACUUM DRYING OVEN

Microprocessor controller

Dual layer tempered glass door for clear observation

Minimum heating time 50% less than traditional vacuum oven

Interior chamber made from stainless steel with a mirror polished finish making for easy cleaning

Standard with vacuum pump

Provides programmability for 7 periods and 9 steps for each period making for a total of 63 programmable step



SPECIFICATIONS

Model	OVE078-913
Capacity	913 L
Electrical requirement	AC3 80 V 50 HZ
Controller	LCD Programmability controller
Power Consumption	5600 W
Temperature range	RT+ 10 - 200°C
Display Resolution	0.1
Temperature Stability	±1°C
Vacuum Degree	133 Pa
Ambient Temperature	+5 - 40°C
Interior Dimension	750Wx1160Dx1050H mm
Exterior Dimension	1400Wx1395Dx2010H mm
Shelves	5 pcs
Chamber Material	304 Stainless steel

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com