

VACUUM DRYING OVEN OVE079-032



VACUUM DRYING OVEN OVE079-032

OVE079-032 VACUUM DRYING OVEN

- Programmable vacuum cycle
- LCD vacuum degree window
- Suitable for drying thermo-sensitive material
- Shorter drying time for sample that does not dry easily
- Dual layer tempered glass door for clear observation
- Fully automatic electromagnetism controlled vacuum
- Minimum heating time 50% less than traditional vacuum oven



SPECIFICATIONS

| | |
|------------------------|---|
| Model | OVE079-032 |
| Capacity | 32 L |
| Electrical requirement | AC 220 50 HZ |
| Temperature Range | RT+10 - 200°C |
| Display Resolution | 0.1°C |
| Temperature Stability | ±1°C |
| Vacuum Gauge | Digital Display |
| Vacuum Degree | 133 Pa |
| Vacuum Sensor | Resistance silicon tube pressure sensor |
| Vacuum Control Range | 10 - 100000 Pa |
| Power Consumption | 1200 W |
| Interior Dimension | 320Wx320Dx320H mm |
| Exterior Dimension | 550Wx490Dx1240H mm |
| Shelves | 2 pcs |
| Chamber Material | 304 Stainless steel |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com