

# VACUUM DRYING OVEN OVE079-090



# VACUUM DRYING OVEN OVE079-090

## OVE079-090 VACUUM DRYING OVEN

Programmable vacuum cycle  
LCD vacuum degree window  
Suitable for drying thermo-sensitive material  
Shorter drying time for sample that does not dry easily  
Dual layer tempered glass door for clear observation  
Fully automatic electromagnetism controlled vacuum  
Minimum heating time 50% less than traditional vacuum oven



### SPECIFICATIONS

Model	OVE079-090
Capacity	91 L
Electrical requirement	AC 220 50 HZ
Temperature Range	RT+10 - 200°C
Display Resolution	0.1°C
Temperature Stability	±1°C
Vacuum Gauge	Digital Display
Vacuum Degree	133 Pa
Vacuum Sensor	Resistance silicon tube pressure sensor
Vacuum Control Range	10 - 100000 Pa
Power Consumption	1350 W
Interior Dimension	450Wx450Dx450H mm
Exterior Dimension	610Wx590Dx1350H mm
Shelves	2 pcs
Chamber Material	304 Stainless steel

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)