

# DRYING OVEN OVE11-040



# DRYING OVEN OVE11-040

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

## OVE11-040 DRYING OVEN

Precision LCD

Cooling fan

Stainless steel inner chamber

Intelligent PID control system

Best quality oven



## SPECIFICATIONS

Model	OVE11-040
Capacity	40 L
Temperature Range	RT+5-250°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±1%
Temperature Uniformity	±1°C
Temperature Resolution	±0.5°C
Timer Range	0-999 min
Internal Dimension	350Wx350Dx320H mm
Exterior Dimension	640Wx550Dx530H mm
Noise Level	≤55 db
Weight	36/39 kg
Power	1000 W
Power Supply	220 V, 50 Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)