

# DRYING OVEN OVE11-070



# DRYING OVEN OVE11-070

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

## OVE11-070 DRYING OVEN

Precision LCD

Cooling fan

Stainless steel inner chamber

Intelligent PID control system

Best quality oven



## SPECIFICATIONS

Model	OVE11-070
Capacity	70 L
Temperature Range	RT+5-250°C
Temperature Accuracy	±0.1°C
Temperature Fluctuation	±1%
Temperature Uniformity	±1°C
Temperature Resolution	±0.5°C
Timer Range	0-999 min
Internal Dimension	450Wx350Dx420H mm
Exterior Dimension	740Wx550Dx630H mm
Package Size	720x740x1470 mm
Weight	47/50 kg
Power	1200 W
Power Supply	220 V, 50 Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)