

DRYING OVEN OVE12-042



DRYING OVEN OVE12-042

Heating Oven, Laboratory Drying Oven

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory general Purpose Oven.

OVE12-042 DRYING OVEN

Micro-computer intelligent control system with LED display

Forced air convection maintains temperature uniformity

High quality stainless steel inner chamber with electrostatic spray treated exterior

Large observation window

Easy to clean semicircular corners with adjustable shelves

Silicon rubber sealing ring prevents heat loss

Over temperature alarm with timer



SPECIFICATIONS

Model	OVE12-042
Capacity	42 L
Temperature Range	RT+5-250°C
Temperature Accuracy	±1°C
Temperature Uniformity	±4.5°C at 250°C
Temperature Resolution	±1°C
Timer Range	0-999 min
Internal Dimension	350Wx350Dx350H mm
Exterior Dimension	530Wx640Dx550H mm
Package Size	650x720x630 mm
Shelves	2
Weight	35/38 kg
Power	1000 W
Power Supply	220V, 50Hz

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com