

DRYING OVEN OVE12-070T



DRYING OVEN OVE12-070T

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

OVE12-070T DRYING OVEN

Digital display all 304 stainless steel inner chamber and outer

Intelligent PID control system

Large observation window

Tight buckle handle



SPECIFICATIONS

Model	OVE12-070T
Capacity	70 L
Temperature Accuracy	±1°C
Temperature Fluctuation	±1%
Time for Heating	≤6 min (20°C to 57°C)
Timer Range	0-999 min
Temperature control accuracy	≤±0.5°C
Internal Dimension	450WX350DX450H mm
Exterior Dimension	740Wx550Dx630H mm
Display Accuracy	0.1°C
Weight	45/48 kg
Power	1200 W
Power Supply	220 V, 50 Hz

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com