

# DRYING OVEN OVE12-225



# DRYING OVEN OVE12-225

Heating Oven, Laboratory Drying Oven

Used in Drying, Roasting, Laboratory, Research, Engineering, Industry, Agricultural, Biology.

Also known as Laboratory Oven, Laboratory general Purpose Oven.

## OVE12-225 DRYING OVEN

Micro-computer intelligent control system with LED display

Forced air convection maintains temperature uniformity

High quality stainless steel inner chamber with electrostatic spray treated exterior

Large observation window

Easy to clean semicircular corners with adjustable shelves

Silicon rubber sealing ring prevents heat loss

Over temperature alarm with timer



## SPECIFICATIONS

Model	OVE12-225
Capacity	225 L
Temperature Range	RT+5-250°C
Temperature Accuracy	±1°C
Temperature Uniformity	±7.0°C at 400°C
Temperature Resolution	±1°C
Timer Range	0-999 min
Internal Dimension	750Wx600Dx500H mm
Exterior Dimension	980Wx965Dx770H mm
Package Size	1100x1045x850 mm
Shelves	2
Weight	93/102kg
Power	2500 W
Power Supply	220V, 60Hz

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)