

# DRYING OVEN OVE12-225T



# DRYING OVEN OVE12-225T

Heating Oven, Laboratory Drying Oven

Used in Drying, Baking, Disinfecting, Sterilization, Laboratory, Research, Engineering, Industry, Thermostat storage, Heating Storage.

Also known as Heating Oven, Laboratory Drying Oven.

## OVE12-225T DRYING OVEN

Digital display all 304 stainless steel inner chamber and outer

Intelligent PID control system

Large observation window

Tight buckle handle



## SPECIFICATIONS

|                         |                         |
|-------------------------|-------------------------|
| Model                   | OVE12-225T              |
| Capacity                | 225 L                   |
| Temperature Accuracy    | $\pm 1^{\circ}\text{C}$ |
| Temperature Fluctuation | $\pm 1\%$               |
| Timer Range             | 0-999 min               |
| Internal Dimension      | 600Wx500Dx750H mm       |
| Exterior Dimension      | 965Wx770Dx980H mm       |
| Weight                  | 94/103 kg               |
| Power                   | 2500 W                  |
| Power Supply            | 220 V, 50 Hz            |

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)